



Universal Ph Test Paper Strip Kit 0-14

Description

- **Universal Ph Test Paper Strip Kit 0-14**
- **Test Range**
- **pH 0-14**
- **Test System**
- **Visual Color Comparison**
- **Application**
- **Testing alkalinity and acid for water, body, saliva and urine**
- **Specification**

- **100 pcs / unit**
- **MOQ**
- **100 units**
- **Shelf Life**
- **New: 2 years from production; Opened: better to use up within 90 days.**
- **Place of Origin**
- **Jilin, China**
- **Instructions**
 - **1. Dip the strip into liquid for 2 seconds and remove it.**
 - **2. Compare against the color chart and get accurate results.**